



Rose Comfort Max tilt and recline chair with casters

HCPCS:









FIGURE A

FIGURE B

FIGURE C

FIGURE D









FIGURE E

FIGURE F

FIGURE G

FIGURE H

- 38 degrees of tilt-lever on the back left positioned easily for the caregiver to adjust
- Recline up to 90 degrees (trendelenburg). Recline lever is Red and located on the back Right
- Seat Width 20" and seat depth 20"
- Seat height adjusts from 17"-22" in 1" increments
- Leg rest: Elevating, length adjustable, lever in the back or side positioned easily for caregiver to adjust
- Footrest: Removable single flip up footrest with angle adjustable foot plate, need also secure so it does not get lost
- Wings: Adjustable Removable for Upper Lateral Support
- Casters: (2) Front 5 total lock (12.5 cm) Rear (1) 5 Directional locking & (1) 5 total lock caster (12.5 cm)
- Standard Padding: Headrest cushion, shoulder bolsters, lower side pads, upper side pads, seat and back pad with cover, calf pad, sole pad

Description
Chair with Caster Full Tilt and Recliner
Description
33"
Blue
Arm Pad: 14" (L) x 3" (W) Leg Pad: 19" (L) x 12.5" (H) Foot Pad: 12" (W) x 22" (H)
23.5"
65"
Up to 90 degrees (trendelenburg)
18"
Standard position: 17" - 22" In full recline: 21.5"
20"
Up to 38 degrees
Lifetime (10 years on frame)
400 lbs
Between Shoulder Width: 22.5"
16 guage tubular steel frame, 1 1/2" vinyl straps, flat head screws, foam
Description
Huntington's Padding Package CLN400-20, EA