Clinical Care Recliner

- The Series FR577RG and FR577RGH Recliners feature three
 positions for patient comfort. The fourth position is a
 Trendelenburg (or Shock) Position that must be activated by a
 caregiver the Trendelenburg Position CANNOT be activated
 by a patient while in the recliner. In addition, the Series FR577RG
 features a back that can be repositioned via the Adjustment
 Lever located on the Push Handle on the back of the recliner.
 - °The first position is the **Upright Position**. This position is ideal for patient transport and can be used to initiate patient treatment / therapy.
 - °The second position is the **First Recline Position** and is also known as the "TV" position. This position is ideal for patient treatment, rest and reading.
 - [°]The third position is the **Second Recline Position** and is also known as the "Full Recline" position. This position is ideal for patient treatment, rest and reading.
- Improved wider and softer armrest for patient comfort.
- Top section of headrest can be moved forward for enhanced patient comfort.

- Softer seat and back upholstery for improved patient comfort.
 Varied density seat back and redesigned seating position provide excellent comfort and support during long term sitting.
- Easy to activate gas cylinder mechanism for quick Trendelenburg activation and improved lift.
- Durable, high-impact side tables sit flush with recliner arm to conveniently hold prep materials. Easy to actuate by caregiver. Side table surface is smooth and features tapered edges for easy cleaning.
- Easily removed side panels enable easy access to recliner for cleaning and maintenance.
- Four Tente® swivel casters facilitate transport and maneuverability in confined areas. Three total locking casters with easy to locate black locking tab lock wheel and swivel with single motion. Directional wheel lock with high visibility green tab allows recliner to track straighter with less effort during transport.
- 290 lb maximum weight capacity, EVENLY DISTRIBUTED.
- Available With Heat and Massage Option.



The New Back Profile, Featuring Specially-Contoured High-Density Foam, Enhances Patient Comfort and Support



Articulating Headrest in Forward Position

Specifications

Overall Height	47"
Overall Width	31"
Overall Depth	35"
Width Between Arms	22"
Arm Height From Floor	25.5"
Arm Height From Seat	8"
Seat Height	22"
Seat Width	21.5"

Seat Depth21"
Depth Fully Reclined
In Trendelenburg74"
Width With Both Side Tables
in Upright Position51"
Shipping Weight 127 lb
Max. Weight Capacity 290 lb

















Colors may vary slightly from actual vinyl

Standard	with Heat & Massage		Standard	with Heat & Massage	
FR577RG401	FR577RGH401	Vintage Gold	FR577RG8615	FR577RGH8615	Berry
FR577RG409	FR577RGH409	Warm Taupe	FR577RG8624	FR577RGH8624	Granny Smith
FR577RG427	FR577RGH427	Blue Ridge	FR577RG8626	FR577RGH8626	Eucalyptus
FR577RG432	FR577RGH432	Imperial Blue	FR577RG8628	FR577RGH8628	Sea
FR577RG454	FR577RGH454	Royal Blue	FR577RG8801	FR577RGH8801	Marine Blue
FR577RG6705	FR577RGH6705	Dolce Sand	FR577RG8805	FR577RGH8805	Jet
FR577RG6714	FR577RGH6714	Dolce Burgundy	FR577RG8808	FR577RGH8808	Champagne
FR577RG6717	FR577RGH6717	Dolce Moss	FR577RG8809	FR577RGH8809	Taupe
FR577RG6724	FR577RGH6724	Dolce Jet	FR577RG8818	FR577RGH8818	Sage
FR577RG6725	FR577RGH6725	Dolce Indigo	FR577RG9201	FR577RGH9201	Cobblestone
FR577RG6727	FR577RGH6727	Dolce Meteor	FR577RG9202	FR577RGH9202	Periwinkle
FR577RG8504	FR577RGH8504	Sand	FR577RG9205	FR577RGH9205	Gypsum
FR577RG8516	FR577RGH8516	Chestnut	FR577RG9207	FR577RGH9207	Ice Blue
FR577RG8533	FR577RGH8533	Port	FR577RG9208	FR577RGH9208	Wineberry
FR577RG8561	FR577RGH8561	Mulberry	FR577RG9210	FR577RGH9210	Cypress
FR577RG8567	FR577RGH8567	Dove	FR577RG9212	FR577RGH9212	Blue Jade
FR577RG8584	FR577RGH8584	Taupe	FR577RG9214	FR577RGH9214	Steel Blue
FR577RG8602	FR577RGH8602	Oatmeal	FR577RG9715	FR577RGH9715	Eucalyptus
FR577RG8604	FR577RGH8604	Concrete	FR577RG9722	FR577RGH9722	Desert Sage
FR577RG8611	FR577RGH8611	Terra Cotta	FR577RG9764	FR577RGH9764	Steel Grey

